




Introduction to Next-Generation Sequencing




Junfan Huang

MRC Cancer Unit
University of Cambridge


CRUK Bioinformatics Summer School 2021
21th July 2021




DNA sequencing



DNA sequencing



Cost per genome




Illumina sequencers



	MiSeq	HiSeq	NovaSeq	Sanger
Reads (millions)	30	3,000	13,000	0.0004
Gigabases/day	7	500	4000	0.001

Illumina sequencers



— = Adenine

— = Thymine

— = Cytosine

— = Guanine

— = Phosphate
backbone


DNA



	MiSeq	HiSeq	NovaSeq	Sanger
Reads (millions)	30	3,000	13,000	0.0004
Gigabases/day	7	500	4000	0.001


Illumina sequencing by synthesis

A. Library Preparation




NGS library is prepared by fragmenting a gDNA sample and ligating specialized adapters to both fragment ends.

B. Cluster Amplification





Library is loaded into a flow cell and the fragments are hybridized to the flow cell surface. Each bound fragment is amplified into a clonal cluster through bridge amplification.

C. Sequencing




Sequencing reagents, including fluorescently labeled nucleotides, are added and the first base is incorporated. The flow cell is imaged and the emission from each cluster is recorded. The emission wavelength and intensity are used to identify the base. This cycle is repeated "n" times to create a read length of "n" bases.





Library Preparation

Add sequencing adapters to DNA fragments





Adapter (20-40 bp)




Flow cell

Library Preparation

Add sequencing adapters to DNA fragments





Cluster generation




Sequencing

Chemistry for sequencing





remove




 terminators

 fluorescent

Sequencing



Chemistry for sequencing





 terminators


 fluorescent

Sequencing



Chemistry for sequencing




 terminators

 fluorescent

Sequencing



Chemistry for sequencing




 terminators


 fluorescent

Sequencing


Chemistry for sequencing






remove




 terminators

 fluorescent

Sequencing





Chemistry for sequencing




 terminators

 fluorescent


Sequencing



A



C





A


...

...


Sequencing




G



T




A



C


Sequencing



Illumina sequencing by synthesis

A. Library Preparation

Genomic DNA 
Fragmentation



NGS library is prepared by fragmenting a gDNA sample and ligating specialized adapters to both fragment ends.

B. Cluster Amplification




Brid



Library is loaded into a flow cell and the fragments are hybridized to the flow cell surface. Each bound fragment is amplified into a clonal cluster through bridge amplification.

C. Sequencing




The emission wavelength and intensity are used to identify the base. This cycle is repeated "n" times to create a read length of "n" bases.



FASTQ

Multiplexing

- Sequences **multiple samples** at **the same time**
- Blocks against possible technical bias caused by differences between flow cell lanes
- Sequences **small genomes** or **specific genomic regions**.



Unaligned sequence: FASTQ

FASTQ header decoded (Illumina example):

	Machine ID	Run	Flow cell ID	Lane	Tile	Tile coordinates		Read		Barcode
						X	Y	Idx	Filter	
Label	@K00359:71:HJJL7BBXX:3:1101:1996:1508 1:N:0:ATCACG									
Sequence	AAAATTCCAAGCTGGTTTCAACAGTACTTTGTTTCCAGAACAAAGAAATG									
Label	+									
Quality score	AAAFFJJJJJJFJJ<J<FJJJJJJJJJJJJJJJJJJFJJFJJJJJFFJFJJJJJJJ<									
	ASCII code									

Worst quality

Best quality

!"#\$%&'()*+,-./0123456789:;<=>?@ABCDEFGHIJKLMNOPQRSTUVWXYZ[\]^_`abcdefghijklmnopqrstuvwxyz{|}~

Alignment



Trimmed DNA sequences

?



GRCh38


Human genome

Alignment



?

GRCh38





10 min break!